Acryl Ort SC

CE 0476

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GB **DIRECTIONS:** Self-curing acrylic resin for orthodontic apparatus, particularly indicated for spray technique.

PREPARATION: It is advisable to work at a room temperature of 22°/23°C.

a) For the pug technique: Mix three parts of powder (polymer) and one part of liquid (monomer). Mix vigorously for about 20 seconds and let it rest for at least 2 minutes. The working time is about 5-6 minutes. After this time, immerse the apparatus in a pressure cooker or a polymerizing machine at 1.5-2 bar for 10 minutes at a temperature of 40°/45°C. If the product isn't subjected to pressurized polymerization, an exothermic reaction will begin after 9 minutes until total polymerization in a free environment.

b) For the direct spray technique: Put the plaster model in water until it is completely saturated to avoid porosity in the finished product. Dry with air and isolate. After having prepared the model to be treated (hooks and arch-scrows with hard wax fixture), sprinkle the polymer directly on the model with a suitable container. Pour a sufficient quantity of the monomer to achieve the completely saturated product. Repeat the operation again until the prosthesis is complete. At the end of the modelling of the apparatus with the resin, it is advisable to remove the excess with a cutting spatula. The working time at this stage is about 5 minutes. Immerse the apparatus in a pressure cooker or polymerizing machine at 1.5-2 bar for 10 minutes at a temperature of 40°/45°C. If the product isn't subjected to pressurized polymerization, an exothermic reaction will begin after 9 minutes until total polymerization in a free environment.

N.B.: If the above procedures are followed, a prosthesis will be achieved with a high stability and resistance to ultimate tensile stress.

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Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction. May cause respiratory irritation. Keep away from heat/sparks/open flames/hot surfaces. - No smoking. Avoid breathing dust/fume/gas/mist/vapours/Smoke. Wear effective gloves/protective clothing/eye protection. Rinse skin with water/shower. Dispose of contents/container in accordance with local/regional/national/international regulations. Product contains: methyl methacrylate, ethylene dimethacrylate.

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